

Image	AcquiredDate	The acquisition date of the Image.
	Description	A multi-line description for the image.
	Name	A short description for the image. This would be used to, for example, select the image from a list.
	ImagingEnvironment	
	AirPressure	AirPressure in millibars[mbar].
	CO2Percent	%CO2 as a percent-fractions from 0.0 to 1.0 [%].
	Humidity	Humidity as a percent-fraction from 0.0 to 1.0 [%].
	Temperature	Temperature [degrees Celsius].
	ObjectiveSettings	
	CorrectionCollar	An adjustable ring on the objective that corrects for changes in immersion medium refractive imdex. Arbitrary scale and unitless.
Medium	A description of a Medium used for the lens. e.g., Oil, Water, WaterDipping, Air, Multi, Glycerol, Other	
RefractiveIndex	Refractive index is that of the immersion medium.	
Pixels		Defines the location and paramater sof the Pixels, the actual binary image data
DimensionOrder	The order in which the individual planes of data are interleaved. e.g., XYZCT, XYZTC, XYCTZ, XYCZT, XYTCZ, XYTZC	
PhysicalSizeX	} Physical size in x, y, z of a pixel in microns[um]	
PhysicalSizeY		
PhysicalSizeZ		
SizeC		
SizeT	} Dimensional size x, y, z, c, t of pixel data array	
SizeX		
SizeY		
SizeZ		
TimeIncrement	Used for time series that have a global timing specification instead of per-timepoint timing info, e.g., a video stream. [s].	
Type	The variable type used to represent each pixel in the image. e.g., int8, int16, int32, uint8, uint16, uint32, float, bit, double, complex, double-complex	
BinData		If the pixel data is stored dierctly in the XML it is enclosed in BinData Elements
TiffData		If the pixel data is stored in an OME-TIFF file it is described by TiffData Elements
FirstC	} The TiffData element describes how the TIFF IFD numbers are mapped to the Pixels.	
FirstT		
FirstZ		
IFD		
PlaneCount		
UUID		The TiffData UUID and Filename are used for multi-file datasets to maintain connections between the files
FileName		
Channel		
AcquisitionMode	AcquisitionMode describes the type of microscopy performed. e.g., WideField, LaserScanningMicroscopy, LaserScanningConfocal, SpinningDiskConfocal, SlitScanConfocal, MultiPhotonMicroscopy, StructuredIllumination, SingleMoleculeImaging, TotalInternalReflection, FluorescenceLifetime, Spectrallmaging, FluorescenceCorrelationSpectroscopy, NearFieldScanningOpticalMicroscopy, SecondHarmonicGenerationImaging, Other	
Color	A color used render this channel	
ContrastMethod	The technique used to achieve contrast. e.g., Brightfield, Phase, DIC, HoffmanModulation, ObliqueIllumination, PolarizedLight, Darkfield, Fluorescence, Other	
EmissionWavelength	Emission wavelength of excitation for a particular channel, in nanometres [nm].	
ExcitationWavelength	Excitation wavelength of excitation for a particular channel, in nanometres [nm].	
Fluor	The name of the fluorophore used to produce this channel.	
IlluminationType	The method of illumination used to capture the channel. e.g., Transmitted, Epifluorescence, Oblique, NonLinear, Other	
Name	A short name for the channel, used to, for example, identify the channel from a list.	
NDFilter	Specifies the combined effect of any neutral density filters used. [% Transmittance]	
PinholeSize	Specifies adjustable pin hole diameters for confocal microscopes (microns [um]).	
PockelCellSetting	Amount the polarization of the beam introduced by Pockel Cell, if any.	
SamplesPerPixel	The number of samples the detector takes to form each pixel value.	
DetectorSettings		
Binning	Represents the number of pixels that are combined to form larger pixels. e.g., 1x1, 2x2, 4x4, 8x8, Other	
Gain	The Gain of the detector.	
Offset	The Offset of the detector.	
ReadOutRate	Detector read speed (MHz)	
Voltage	The Voltage of the detector. volts[V]	
LightSourceSettings		
Attenuation	The Attenuation of the light source. [%]	
Wavelength	The Wavelength of the light source. [nm]	
Plane		
DeltaT	Elapsed time since the beginning of the experiment [s]	
ExposureTime	Elapsed time during image recording. [s]	
PositionX	} The x, y, z position of the stage. [µm]	
PositionY		
PositionZ		
StageLabel		
Name	Short name for this stage location This would be used to, for example, identify the channel from a list.	
X	} The labeled x, y, z position of the stage. [µm]	
Y		
Z		