Install - Web

OMERO.web Install and Configuration

OMERO.web is the web application component of the OMERO platform which allows for the management, viewing (a fully multi-dimensional image viewer), attachment and association of files to images, datasets, and projects from one or more OMERO.server's. It also includes webadmin for managing users and groups. With one or two commands, OMERO.web can be started with the Django lightweight development Web server bound to port 8000 on 127.0.0.1 when you deploy your OMERO.server instance.

We **do not** recommend using this server in a production setting. It has not gone through security audits or performance tests. Apache with mod_python is the preferred setup for deploying a Django application on a production server. You can find more information on the Django site here.

Prerequisites

- OMERO 4 and its prerequisites
- If you are using Python 2.4, you will need pysqlite. Use version 2.0.3 or higher. Python 2.5 ships with an SQLite wrapper in the standard library, so you do not need to install anything extra in that case
- Python Imaging Library and Matplot Lib
 - PIL and Matplot Lib are an optional requirement for drive space utilization graphs
- Apache with mod python

PYTHONPATH

on non-Windows platform

```
export PYTHONPATH=$PYTHONPATH:/home/omero/omero dist/lib/python/
```

```
C:\>set PYTHONPATH=%PYTHONPATH%;C:\OMERO\lib\python
C:\>echo %PYTHONPATH%
C:\Ice-3.3.0\python;C:\omero-Beta4.0.0\lib\python
```

Configuration for Django Web server (non-production)

On non-Windows* platforms, OMERO.web is configured to be started with a single com

Start up OMERO.web manually:

```
$ cd /home/omero/omero-dist/
$ bin/omero admin ice
Ice 3.3.0 Copyright 2003-2008 ZeroC, Inc.
>>> server enable Web
>>> server start Web
  >>> exit
```

Shut down OMERO.web manually

```
$ cd /home/omero/omero-dist/
$ bin/omero admin ice
Ice 3.3.0 Copyright 2003-2008 ZeroC, Inc.
>>> server disable Web
>>> server stop Web
>>> exit
```

■ Or have it started automatically

If you wish to have the server always on, configure the following in /home/omero/omero_dist/etc/grid/default.xml:

```
<server-instance template="ShellTemplate" id="Web" act="always"/>
```

■ If you require to configure the out of the box setup to listen for webadmin and webclient connections on different host configure the following:

Configure the following in /home/omero/omero_dist/etc/grid/templates.xml by adding after --norcload:

```
</server-template>
```

In your /home/omero/omero dist/lib/python/omero/plugins/directory configure the following in server.py:

```
django = ["python", "manage.py", "runserver", "0.0.0.0:8080", "--noreload"]+list(args)
```

On Windows* platforms, please follow the instruction on OMERO install - windows

Configuration Guide for mod_python (production)

■ If you've enabled the OMERO.web instance that is using the Django lightweight development web server, you will need to shut it down:

Shut down the OMERO.web manually

```
$ cd /home/omero/omero_dist/
$ bin/omero admin ice
Ice 3.3.0 Copyright 2003-2008 ZeroC, Inc.
>>> server disable Web
>>> server stop Web
```

In your /home/omero/omero dist/etc/qrid/ directory re-configure the Web server in default.xml to be on-den

1 of 3 21/10/2009 12:34

```
<server-instance template="ShellTemplate" id="Web" act="on-demand"/>
```

■ In your /home/omero/omero_dist/lib/python/omeroweb/ directory configure the following in settings.py:

■ DATABASE NAME

```
DATABASE_NAME = '/home/omero/webdb/db.sqlite3'
```

■ LOGDIR

```
# LOGS
# to change the log place, please specify new path for LOGDIR.
# to change the log place, please specify new path for LOGDIR.
# LOGDIR = os.path.join(os.path.dirname(_file_), 'log')
# LOGDIR = os.path.join(os.path.join(os.path.join(os.path.dirname(_file_), '...'), 'var'), 'log')
LOGDIR = '/homme/ommero/weblog/'
```

■ Notification options: EMAIL_NOTIFICATION, EMAIL_SENDER_ADDRESS and EMAIL_SMTP_SERVER

```
# Notification
# Application allows to notify user about new shares
EMAIL NOTIFICATION = True
EMAIL SONDER ADDRESS = 'omero-admin@example.com'
EMAIL_SMTP_SERVER = 'smtp.example.com'
# Optional settings
# EMAIL_SMTP_PORT = 25
# EMAIL_SMTP_USER = 'login'
# EMAIL_SMTP_USER = 'login'
# EMAIL_SMTP_TLSE = True
```

· Notification option for sharing:

```
# If web server is running behind a proxy server please set the host.
# That option is required by share notification sendere.
APPLICATION_HOST='http://www.domain.com:80'
```

■ Create log and database directories with apache user:apache group ownership

```
mkdir /home/omero/weblog
chown apache_user:apache_group /home/omero/weblog
mkdir /home/omero/webdb
chown apache_user:apache_group /home/omero/webdb
```

■ Initialize the Django settings database, specifying an OMERO.web superuser username, password and e-mail address:

```
# cd /home/omero/omero_dist/lib/python/omeroweb/
# sudo -u apache_user PYTHONPATH=/home/omero/omero_dist/lib/python/ python manage.py syncdb

Creating table django_admin_log
Creating table auth_permission

Creating table auth_user
Creating table auth_user
Creating table auth_message
Creating table auth_message
Creating table django_content_type
Creating table django_session
Creating table django_siste
Creating table webadmin_gateway

You just installed Django's auth system, which means you don't have any superusers defined.
Would you like to create one now? (yes/no): yes
Username (Leave blank to use 'omero': omero
E-mail address: omero-admin@example.com
Password:
Password (again):
Superuser created successfully.
Installing custom SQL for webadmin.Gateway model
Installing index for auth.Permission model
Installing index for auth.Permission model
Installing index for auth.Message model
```

■ Configure mod python:

```
CLocation "/">
SetHandler python-program
PythonHandler django.core.handlers.modpython
SetEnv DJANGO_SETTINGS_MODULE omeroweb.settings
PythonDebug On

(for 4.0.0) PythonPath "['/home/omero/omero_dist/lib', '/home/omero/omero_dist/lib/python', '/home/omero/omero_dist/lib/omeroweb'] + sys.path"

(as of 4.0.1) PythonPath "['/home/omero/omero_dist/lib/python', '/home/omero/omero_dist/lib/python/omeroweb'] + sys.path"

</Location>
```

- Go to the link http://your_host/webclient/ or http://your_host/webadmin/
- If you wish to modify the predefined list of the OMERO.servers please go to admin panel: http://your_host/admin/webadmin/gateway/. To login, use the user superuser username and password created you created when initializing the Django settings database above.

Once you have deployed and started the server you can use your browser to access the OMERO.webadmin or OMERO.webclients

- http://your_host/webadmin
- http://your_host/webclient



2 of 3 21/10/2009 12:34



Troubleshooting

My OMERO install doesn't work! What do I do now!?! Examine the Troubleshooting page and if all else fails post a message to our ome-users mailing list discussed on the Community page.

3 of 3 21/10/2009 12:34